

Indi-Strecker




Features

Sizes available for the **smallest to large** patients

Changing position of the ball joint adjusts flexion angle

+ **Force unit** assists knee extension and helps when sitting down and standing up



Standard


Fully assembled complete with everything needed for a functional add-on to the Indi knee joint

Accessory available as **combi-product** with the Indi knee joint


Serial number for traceability

Clearly illustrated instructions for use


Fig. left




We recommend using our configurator to make your selection.
The configurator considers, among other factors, weight to height ratios.



Weight Classification (body weight)

Article No. Combi-product		KAFO	HKAFO		Article No.
Indi-Strecker-22		not yet available	not yet available		Indi-Strecker--22
Indi-Strecker-25-Plus		not yet available	not yet available		Indi-Strecker-25-Plus
Indi-Strecker-28		not yet available	not yet available		Indi-Strecker-28
Indi-Strecker-30		up to 50 kg	up to 50 kg		Indi-Strecker-30
Indi-Strecker-30-Plus		up to 70 kg	up to 70 kg		Indi-Strecker-30-Plus
Indi-Strecker-34		not yet available	not yet available		Indi-Strecker-34
Indi-Strecker-34-Plus		not yet available	not yet available		Indi-Strecker-34-Plus

								
Locked	Extension 5°	Flexion 115°	Posterior Offset	Modular	Force Options	No Grinding Necessary	Noise Dampening	Water Resistant