

Klassik-Konso Slide



Features

- + Wear-free, **magnetic** safety release
- + Ergonomic **plastic grip**
- + **Pre-select system** opens only when knee is in full extension and drop lock is relieved

Cable with or without housing




Fig. left


Standard

Fully assembled as one item for clear cost estimates and uncomplicated ordering


Accessory available as **combi-product** with the Klassik-Konso knee joint

Serial number for traceability


Clearly illustrated instructions for use



We recommend using our configurator to make your selection.
The configurator considers, among other factors, weight to height ratios.




Weight Classification (body weight)


Article No. Combi-product		KAFO		Article No.
Klassik-Konso-Slide-30		up to 35 kg		Slide-N-30
Klassik-Konso-Slide-32		up to 50 kg		Slide-N-32
Klassik-Konso-Slide-36		up to 70 kg		Slide-N-36
Klassik-Konso-Slide-40		up to 85 kg		Slide-N-40
Klassik-Konso-Slide-40-Plus		up to 85 kg		Slide-N-40-Plus
Only available as replacement item				



Locked



Extension



Flexion

Posterior Offset

Modular

Force Options



No Grinding Necessary

Noise Dampening

Water Resistant