

Konso-Move-S



Features




Fig. left


Standard

Fully assembled complete with everything needed for a functional add-on to the Konso knee joint


Available **separately or as combi-product** together with a Konso knee joint

Serial number for traceability

Clearly illustrated instructions for use



We recommend using our configurator to make your selection.
The configurator considers, among other factors, weight to height ratios.



Weight Classification (body weight)

Article No. Combi-product		KAFO		Article No.
Konso-Move-S-22		up to 20 kg		Move-S-22
Konso-Move-S-25		up to 30 kg		Move-S-25
Konso-Move-S-28		up to 35 kg		Move-S-28
Konso-Move-S-30		up to 50 kg		Move-S-30
Konso-Move-S-32		up to 70 kg		Move-S-32
Konso-Move-S-36		up to 85 kg		Move-S-36
Konso-Move-S-40		up to 85 kg		Move-S-40

Konso contoured | Konso straight

								
Locked	Extension	Flexion	Posterior Offset	Modular	Force Options	No Grinding Necessary	Noise Dampening	Water Resistant