

Stehfeld Standing Platform



Features

Sizes available for the **smallest to large** patient

Unrivalled **20-month guarantee**

Engraved **serial number** for traceability

Clearly illustrated instructions for use

Standard

Fully assembled as one item for clear cost estimates and uncomplicated ordering

+ Tilting connecting plate for **posture control exercises** during therapy

+ **Silicone bumper** to avoid damages and protect against injury


+ **Wheel locks** keep podium stationary


+ **Pyramid adapters** for precise adjustment

+ **Height adjustable** to accommodate growth

+ Unscrew **2 grips** to detach body shell and posture plate from the base


+ Applies the Helga Hirschfeld therapeutic principle to **induce postural reactions** by displacing the body and prompting a return to equilibrium






Article No.	Height Classification	Base (l x w)	Posture Plate	Foot Length	Foot to Foot (lateral)
Stehfeld-klein	60 cm to 100 cm	40 cm x 40 cm	22,5 cm x 29 cm	12,5 cm to 17 cm	17,0 cm to 23 cm
Stehfeld-mittel	80 cm to 130 cm	45 cm x 45 cm	25,5 cm x 36 cm	14,5 cm to 22 cm	20,5 cm to 31 cm
Stehfeld-groß	120 cm to 160 cm	50 cm x 50 cm	28,5 cm x 42 cm	16,5 cm to 25 cm	22,5 cm to 35 cm
Stehfeld-groß-breit	120 cm to 160 cm	50 cm x 70 cm	28,5 cm x 42 cm	16,5 cm to 25 cm	22,5 cm to 35 cm

With anti-tip protection to prevent toppling backwards




Article No. Combi-product	Height Classification	Base (l x w)	Posture Plate	Foot Length	Foot to Foot (lateral)
Stehfeld- KS -klein	60 cm to 100 cm	40 cm x 40 cm	22,5 cm x 29 cm	12,5 cm to 17 cm	17,0 cm to 23 cm
Stehfeld- KS -mittel	80 cm to 130 cm	45 cm x 45 cm	25,5 cm x 36 cm	14,5 cm to 22 cm	20,5 cm to 31 cm
Stehfeld- KS -groß	120 cm to 160 cm	50 cm x 50 cm	28,5 cm x 42 cm	16,5 cm to 25 cm	22,5 cm to 35 cm
Stehfeld- KS -groß-breit	120 cm to 160 cm	50 cm x 70 cm	28,5 cm x 42 cm	16,5 cm to 25 cm	22,5 cm to 35 cm

With anti-tip pro



Article No. Combi-product	Height Classification	Base (l x w)	Posture Plate	Foot Length	Foot to Foot (lateral)
Stehfeld- HG -klein	60 cm to 100 cm	40 cm x 40 cm	22,5 cm x 29 cm	12,5 cm to 17 cm	17,0 cm to 23 cm
Stehfeld- HG -mittel	80 cm to 130 cm	45 cm x 45 cm	25,5 cm x 36 cm	14,5 cm to 22 cm	20,5 cm to 31 cm
Stehfeld- HG -groß	120 cm to 160 cm	50 cm x 50 cm	28,5 cm x 42 cm	16,5 cm to 25 cm	22,5 cm to 35 cm
Stehfeld- HG -groß-breit	120 cm to 160 cm	50 cm x 70 cm	28,5 cm x 42 cm	16,5 cm to 25 cm	22,5 cm to 35 cm

With hand grips and anti-tip protection




Article No. Combi-product	Height Classification	Base (l x w)	Posture Plate	Foot Length	Foot to Foot (lateral)
Stehfeld- HG-KS -klein	60 cm to 100 cm	40 cm x 40 cm	22,5 cm x 29 cm	12,5 cm to 17 cm	17,0 cm to 23 cm
Stehfeld- HG-KS -mittel	80 cm to 130 cm	45 cm x 45 cm	25,5 cm x 36 cm	14,5 cm to 22 cm	20,5 cm to 31 cm
Stehfeld- HG-KS -groß	120 cm to 160 cm	50 cm x 50 cm	28,5 cm x 42 cm	16,5 cm to 25 cm	22,5 cm to 35 cm
Stehfeld- HG-KS -groß-breit	120 cm to 160 cm	50 cm x 70 cm	28,5 cm x 42 cm	16,5 cm to 25 cm	22,5 cm to 35 cm

Body Shell with Posture Plate

Custom body shell used in conjunction with the Stehbrett supports primary deformities and selectively counters any possible secondary deformities.


Only a precise footbed and proper support of the ankle joint enable correct positioning of the knee and hip joints.




Amount of tilt can be individually adapted to the user by positioning supports on the undersurface.




Growable




Stability




Protective Bumper




3-Point Fixation




Removable



Postural Control Exercises



HG = Hand Grips



KS = Anti-tip Protection