













































































Functional Overview Walking Podiums

Podiums are combined with an individual, custom-fabricated body shell to form the walking unit. The body shell is affixed to a universal connecting plate, which is fully interchangeable with all sizes as well as the Stehbrett platform (see "Standing Podiums"). The walking podiums have a low center of gravity to increase stability. Hand grips and anti-tip options are available which can be added at any time and are compatible with all podiums. The Leicht-Walker utilizes carbon parts to reduce weight by approx. 30% compared to the standard aluminum Foot-Walker.

		Individual Adjustments			Body Shell			Modular Accessories			
Noise Dampening	Growable	Vertical Axis	Step Width	Step Length	3-Point Fixation	Removable	Reduced Weight	Hand Grips	Anti-tip Protection		
											Foot-Walker
											Foot-Walker-KS
											Foot-Walker-HG
											Foot-Walker-HG-KS
											Leicht-Walker
											Leicht-Walker-KS
											Leicht-Walker-HG
											Leicht-Walker-HG-KS