





























































Functional Overview Standing Podiums

Podiums are combined with an individual, custom-fabricated body shell to form the standing unit. The body shell is affixed to a universal connecting plate, which is fully interchangeable with all Stehbrett sizes, as well as Foot-Walker and Leicht-Walker podiums (see “Walking Podiums”). The Stehfeld podium comes with an adjustable, tilting connecting plate, which can be used for posture control exercises to augment therapy sessions. Hand grips and anti-tip options are available which can be added at any time and are compatible with all podiums.

			Body Shell			Modular Accessories			
Growable	Stability	Protective Bumper	3-Point Fixation	Removable	Posture Control Exercises	Hand Grips	Anti-tip Protection		
									Stehbrett
									Stehbrett -KS
									Stehbrett -HG
									Stehbrett -HG-KS
									Stehfeld
									Stehfeld -KS
									Stehfeld -HG
									Stehfeld -HG-KS