





























































Functional Overview Ankle Joints

All ankle joints are shipped fully assembled. The joints are free of play, self-lubricating and dirt resistant. Force options are differentiated as follows: blue (light), green (mild), grey (medium) and yellow (strong). **The modular options Turbo, Uno-front, Uno-back, -Flex and -Dor can be added or exchanged on the Shuttle joint at any time.**



Adjustable Alignment	Modular	Force Options	Range of Motion	Swing Phase	Quick Adjustment	Interchangeable Stops	No Grinding Necessary	Noise Dampening	Water Resistant	
										Shuttle
										Shuttle with Turbo
										Shuttle with Uno-front
										Shuttle with Uno-back
										Shuttle with Dor
										Shuttle with Flex
										Easy-T

Adjustable Alignment	Modular	Force Options	Range of Motion	Swing Phase	Quick Adjustment	Interchangeable Stops	No Grinding Necessary	Noise Dampening	Water Resistant		
											Easy-A
											Ela-Dor
											Ela-FI-Ex
											Ela-FI-Ex with Post
											Naveta
											Jumper
											Basic